

Bioject's new **ZetaJet™** system is the latest advance in needle-free delivery systems, offering a significant array of unique product features and patient benefits that can provide highly competitive differentiation to a wide variety of injectable market segments.



Bioject® ZetaJet™

The **Bioject® ZetaJet™** Needle-free Injection System is a compact spring-powered, durable device. It is intended to deliver vaccines and injectable medication either subcutaneously or intramuscularly.

The **Bioject® ZetaJet™** is indicated for both professional use and home use for patients who self-inject.

The auto-disable feature eliminates the risks of re-use of a syringe and prevents possible contamination of syringes and vials.



Bioject® ZetaJet™ Syringe

Studies show a strong preference for a Bioject system versus traditional injection delivery.

**Bioject Inc.**

7180 SW Sandburg St., Ste. 100 Tigard, OR 97223 USA

503-692-8001 [www.bioject.com](http://www.bioject.com)

Customer Service: 800-683-7221 ext. 4122

[customerservice@bioject.com](mailto:customerservice@bioject.com)

## ADVANTAGES

- All three injection depths (SC, IM, ID\*)
- Custom molded exterior, provides a multitude of ergonomic options, to achieve specific attributes for different clinical or branding applications
- Injection volumes range from 0.05 mL to 0.5 mL
- Eliminates needle-phobia thereby enhancing patient experience
- Leverages Bioject's unique, patented pressure profile
- Indicated for both professional and patient self-injection
- Polycarbonate syringe provides exceptional strength, clarity and heat resistance
- Reliable stainless steel inner core enhances durability in rough environments



## CONVENIENCE

- Easier to train and use needle-free systems
- Delivers injection in a fraction of a second
- Reduces sharps waste
- The plunger is pre-assembled into the syringe and can be used for reconstitution and other pre-injection tasks

## SAFETY

- No risk of needle-stick injury
- Elimination of cross-contamination
- Auto-disabling syringe

### Bioject Inc.

7180 SW Sandburg St., Ste. 100 Tigard, OR 97223 USA

503-692-8001 [www.bioject.com](http://www.bioject.com)

Customer Service: 800-683-7221 ext. 4122

[customerservice@bioject.com](mailto:customerservice@bioject.com)

\* Investigational use only